

**Disposable Surgical Mask (Ear Loop and Head Ties) Type IIR**  
**Product Standards: EN14683:2019+AC:2019 ASTM Level 3**




**Head Ties Mask** Issue date 03-06-2020

**Product Type:** Disposable Surgical Mask

**Product Description:** The disposable masks fulfil the requirements of EN 14683:2019+AC:2019 compatible with the medical devices directive 93/42/EEC. Also complies with ASTM Level 3 Type 2R (160mm HG spray resistance).

■ **PRODUCT DESCRIPTION**

Product code	WS-MED2-L3 (EL or HT)
Type	TYPE IIR <input checked="" type="checkbox"/>
Product Standards	EN14683:2019+AC:2019 ASTM Level 3, Type 2R (160mm HG SR)
Colour	White <input type="checkbox"/> Blue <input checked="" type="checkbox"/> Grey <input type="checkbox"/>
Shape	Disposable mask (Ear loop and Head Ties types)
Exhalation Valve	No
Usage	General Protection
Availability	Single-Use Only (2)
Country of Origin	China
Manufacturer / Supplier	 Sichuan Zhengning Medical Equipment Co., Ltd.

**Manufacturer:** Sichuan Zhengning Medical Equipment Co., Ltd.

Address: No. 3, Building 2, No. 388, Food Avenue, Juntun Town, Xindu District, Chengdu, Sichuan, China

**Distributor:** DTS International Group P/L. Australia, USA, China, Singapore, Hong Kong, Africa

Address: Level 1, 414 Victoria Ave., Chatswood, NSW 2818

## Disposable Surgical Mask (Ear Loop and Head Ties) Type IIR Product Standards: EN14683:2019+AC:2019 ASTM Level 3

### ■ MATERIAL

These face masks are free of natural rubber latex and silicon

Outer layer	PP non-woven fabric
Filtration layer	Polypropylene (melt-blown )
Inner layer	PP non-woven fabric
Ear Loop	Polyester
Nose clip	Polyethylene covered steel wire

### ■ PHYSICAL PROPERTIES

Bacterial filtration efficiency	BFE≥98%
Differential pressure	38.2 Pa/cm <sup>2</sup>
Microbial cleanliness	≤30CFU/g
Splash resistance pressure	≥16.0
Materials and construction	Pass

### ■ Weight & Size

Weight	Complete face mask : 4.0 grams
Dimensions	175 x 95 mm
Mask thickness	1.0 mm
Elastic bands dimensions	200×2mm

### ■ Usage

The mask is designed for up to 4 hours of usage. Do not re-use; reuse can result in bacteria growing in the mask; should bacteria begin to develop in the mask, you may notice an unusual odor.

**Manufacturer:** Sichuan Zhengning Medical Equipment Co., Ltd.

Address: No. 3, Building 2, No. 388, Food Avenue, Juntun Town, Xindu District, Chengdu, Sichuan, China

**Distributor:** DTS International Group P/L. Australia, USA, China, Singapore, Hong Kong, Africa

Address: Level 1, 414 Victoria Ave., Chatswood, NSW 2818

**Disposable Surgical Mask (Ear Loop and Head Ties) Type IIR**  
**Product Standards: EN14683:2019+AC:2019 ASTM Level 3**

**GTTC**



**Test Report**

(Electronic version)

**No: 20R000277**



**Guangzhou Inspection Testing and Certification Group Co., Ltd.**  
Add: No. 1, Zhujiang Road, Panyu District, Guangzhou, Guangdong, P.R.China.

Tel: +86-20-61994598 61994599

**Manufacturer:** Sichuan Zhengning Medical Equipment Co., Ltd.

Address: No. 3, Building 2, No. 388, Food Avenue, Juntun Town, Xindu District, Chengdu, Sichuan, China

**Distributor:** DTS International Group P/L. Australia, USA, China, Singapore, Hong Kong, Africa

Address: Level 1, 414 Victoria Ave., Chatswood, NSW 2818

**Disposable Surgical Mask (Ear Loop and Head Ties) Type IIR**  
**Product Standards: EN14683:2019+AC:2019 ASTM Level 3**

■ **Packaging (Box version)**

Box	Size: 200x100×88 mm    Weight: 40g
Carton	Size: 460x520×420 mm    Weight: 1.3kg
Packaging	50pcs/Box    40Box/CTN
Total	2000pcs    Gross Weight: 13.3kg
Shelf Life	2 years
Storage Instructions	Store in a ventilated, dark and dry place, away from fire and pollution



**Manufacturer:** Sichuan Zhengning Medical Equipment Co., Ltd.

Address: No. 3, Building 2, No. 388, Food Avenue, Juntun Town, Xindu District, Chengdu, Sichuan, China

**Distributor:** DTS International Group P/L. Australia, USA, China, Singapore, Hong Kong, Africa

Address: Level 1, 414 Victoria Ave., Chatswood, NSW 2818